

PROTECT YOUR STAFF FROM THE HEALTH DANGERS INVOLVED FROM TAKING RADIOGRAPHS

Why Monitor Exposure to X-ray Radiation

- OSHA and many state agencies require the use of x-ray monitoring devices for everyone involved in the use of any radiation-emitting equipment
 - Failure to comply with these laws could result in fines
- Historical data from x-ray monitoring devices can be invaluable in disproving liability lawsuits regarding alleged exposure

X-ray Monitoring

- The PL Medical x-ray badge service uses state-of-the-art dosimetry badges. Our badges detect exposure with four elements and report exposure at three tissue depths (deep, shallow & eye)
- The lightweight badges are approximately 1" x 2" and clip easily onto your clothing

Quarterly Badge Service

- Each person in your clinic enrolled in the x-ray monitoring service will automatically receive a new badge every 3 months
- Reply envelopes are included, allowing you to quickly return the old badges for processing—just add postage and ship via USPS
- A printed report will be mailed to you for your records and are available online for immediate viewing

Service and Quality You Expect

- If any reading should ever exceed NRC safety standards you will be contacted by phone immediately
- Access to past reports
- **Lock-In Rates:** PL Medical can guarantee your badge service charges for a 2-year period

No Hidden Fees

- No setup charges
- No charge to modify service due to staff changes (immediate and full credit for individuals removed from the service)
- No charge for replacement copies of reports

If you are currently enrolled in another service, PL Medical can take early enrollments and begin your service when your current service expires. Your account will not be billed until the service through PL Medical begins.

