

- Wording for the top of the page: "PL Medical offers i.v. extension sets with exceptional quality and competitive value. PL Medical I.V. extension sets are distributed in over 50 countries worldwide. Our manufacturing facilities are both ISO 9001 and ISO 13485 certified. Our I.V. extension sets have been manufactured with the invaluable inputs from scientists and research and development professionals throughout the world.

Part Number: IT0101

I.V. INFUSION SET, VENTED 73" WITH FLASH BACK BULB


Approximate priming volume – 12 ml

Fluid path sterile

Latex free

20 drops/ml ± 0.1ml

- (1) Protective cover (2) Spike (3) Bacterial retentive air filter (4) Drip chamber (5) Tube (6) Roller clamp (7) Flashback bulb (8) Male adaptor (9) Protective cap



I.V. INFUSION SET, VENTED 78" WITH INJECTION SITE

Part Number: IT0102


Approximate priming volume – 13 ml

Fluid path sterile

Latex free

20 drops/ml ± 0.1ml

- (1) Protective cover (2) Spike (3) Bacterial retentive air filter (4) Drip chamber (5) Slide clamp (6) Roller clamp (7) Tube (8) Injection site (9) Male adaptor (9) Protective cap



I.V. INFUSION SET, MICRODRIP 70" WITH INJECTION SITE AND OPTION LOCK

Part Number: IT0103


Approximate priming volume – 12.2 ml

Fluid path sterile

Latex free

60 drops/ml ± 0.1ml

- (1) Protective cover (2) Spike (3) Filter vent cover (4) Microdrip (5) Drip chamber (6) Slide clamp (7) Roller clamp (8) Tube (9) Injection site (10) Option lock (11) Male adaptor (12) Protective cap



I.V. INFUSION SET, CONVERTIBLE PIN – 78" WITH INJECTION SITE AND OPTION LOCK

Part Number: IT0104


Approximate priming volume – 12 ml

Fluid path sterile

Latex free

20 drops/ml \pm 0.1ml

(1) Protective cover (2) Spike (3) Filter vent cover (4) Drip chamber (5) Slide clamp (6) Roller clamp (7) Tube (8) Injection site (9) Option lock (10) Male adaptor (11) Protective cap



I.V. INFUSION SET, CONVERTIBLE PIN – 106” WITH BACK CHECK VALVE AND 3 INJECTION SITE AND OPTION LOCK

Part Number: IT0106


Approximate priming volume – 17 ml

Fluid path sterile

Latex free

20 drops/ml \pm 0.1ml

(1) Protective cover (2) Spike (3) Filter vent cover (4) Drip chamber (5) Backcheck valve (6) Injection site (7) Slide clamp (8) Roller Clamp (9) Tube (10) Option lock (11) Male adaptor (12) Protective cap



I.V. INFUSION SET WITH DUO VENT SPIKE 103” INJECTION SITE MALE LUER LOCK ADAPTOR RETRACTABLE COLLAR

Part Number: IT0107


Approximate priming volume – 17 ml

Fluid path sterile

Latex free

20 drops/ml \pm 0.1ml

(1) Protective cover (2) Spike (3) Filter vent cover (4) Drip chamber (5) Slide clamp (6) Roller clamp (7) Tube (8) Injection site (9) Option lock (10) Male adaptor (11) Luer cap



I.V. INFUSION SET WITH DUO VENT SPIKE 104” 2 INJECTION SITES MALE LUER LOCK ADAPTOR RETRACTABLE COLLAR

Part Number: IT0108


Approximate priming volume – 17 ml

Fluid path sterile

Latex free

20 drops/ml \pm 0.1ml

(1) Protective cover (2) Spike (3) Filter vent cover (4) Drip chamber (5) Slide clamp (6) Roller clamp (7) Tube (8) Injection site (9) Option lock (10) Male adaptor (11) Protective cap



I.V. INFUSION SET 78", VENTED WITH 3 INJECTION SITES


Part Number: IT0114

Approximate priming volume – 13 ml

Fluid path sterile

Latex free

(1) Protective cover (2) Spike (3) Bacterial retentive air filter (4) Drip chamber (5) Check valve (6) Tube (7) Cap (8) Injection site (9) Slide clamp (10) Roller clamp (11) Male adaptor (12) Luer cap



I.V. INFUSION SET 78", VENTED WITH ONE INJECTION SITE


Part Number: IT0117

Approximate priming volume – 13 ml

Fluid path sterile

Latex free

(1) Protective cover (2) Spike (3) Bacterial retentive air filter (4) Drip chamber (5) Slide clamp (6) Roller clamp (7) Tube (8) Cap (9) Injection site (10) Male adaptor (12) Protective cap



I.V. INFUSION SET 78", VENTED WITH 2 INJECTION SITES


Part Number: IT0118

Approximate priming volume – 13 ml

Fluid path sterile

Latex free

(1) Protective cover (2) Spike (3) Bacterial retentive air filter (4) Drip chamber (5) Slide clamp (6) Tube (7) Roller clamp (8) Injection site (9) Male adaptor (10) Protective cap (11) Cap



I.V. INFUSION SET 78", VENTED WITH 3 INJECTION SITES


Part Number: IT0119

Approximate priming volume – 13 ml

Fluid path sterile

Latex free

(1) Protective cover (2) Spike (3) Bacterial retentive air filter (4) Drip chamber (5) Tube (6) Cap (7) Injection site (8) Slide clamp (9) Roller clamp (10) Male adaptor (11) Protective cap



I.V. INFUSION SET 78", VENTED WITH 0.2 MICRON I.V. LINE FILTER AND 2 INJECTION SITES


Part Number: IT0120

Approximate priming volume – 13 ml

Fluid path sterile

Latex free

(1) Protective cover (2) Spike (3) Bacterial retentive air filter (4) Drip chamber (5) Tube (6) Cap (7) Injection site (8) Slide clamp (9) Roller clamp (10) I.V. filter 0.2 micron (11) Male adaptor (12) Protective cap



I.V. INFUSION SET 78", VENTED WITH CHECK VALVE AND TWO INJECTION SITES

Part Number: IT0121

Approximate priming volume – 13 ml

Fluid path sterile

Latex free

(1) Protective cover (2) Spike (3) Bacterial retentive air filter (4) Drip chamber (5) Check valve (6) Tube (7) Cap (8) Injection site (9) Slide clamp (10) Roller clamp (11) Male adaptor (12) Protective cap

