

The image features two different types of IV administration sets against a purple background with a white grid pattern. On the left is a standard IV set with a Y-shaped top, a blue roller clamp, and a blue drip chamber. On the right is a coiled IV set with a large, flexible, coiled section in the middle. Both sets have clear plastic tubing and white connectors at the top and bottom.

IV Line Administration Sets & Extensions

wide range of standard
and custom
configurations available

IV Administration Set

Non-Pyrogenic, Sterile Fluid Path

IT 0109-PL
Injection Port

60 Drops/ml
70in (178cm)

IT 0107-PL
Injection Port

20 Drops/ml
104in (264cm)

IT 0114-PL
Three Needleless
Injection Ports &
Check Valve






20 Drops/ml
78in (198cm)

IT 0117-PL
Needleless
Injection Port

20 Drops/ml
78in (198cm)

IT 0118-PL
Two Injection Ports

20 Drops/ml
78in (198cm)



Injection Sets

Latex Free, DEHP Free

Injection Ports

20 Drops/ml
78in (198cm)

IT 0119-PL
Three Injection
Ports

20 Drops/ml
78in (198cm)

IT 0120-PL
Two Needleless
Injection Ports &
0.2 micron Filter





20 Drops/ml
78in (198cm)

IT 0122-PL
Two Needleless
Injection Ports

20 Drops/ml
78in (198cm)

IT 0360-PL
360in Coiled MRI
Extension Set

20 Drops/ml
78in (198cm)



S1: Spike Assy w/open-close Vent

C1: Roller Clamp

C2: Pinch Clamp

C3: Slide Clamp

S1a: Protective Cap for Spike

J1: Injection Port

J2: Needleless Injection Port

J3: T-Injection Port

S2: Spike w/open-close Vent

L1: Female Luer Lock

L1a: Protective Cap for Female Luer

K1: 3-Way Stopcock

S2a: Protective Cap for Spike

L2: Male Luer Lock

L2a: Protective Cap for Male Luer

M2: Measured Volume Chamber

IV Extension Sets

Non-Pyrogenic, Sterile Fluid Path, Latex Free, DEHP Free

IT 0109-PL
Roller Clamp Line

30in (76cm)

IT 0110-PL
Pinch Clamp Line

30in (76cm)

IT 0111-PL
Two Injection Ports

42-in (106cm)


IT 0112-PL
T-Injection Port

5.25in (13cm)

IT 0113-PL
T-Injection Port

5.25in (13cm)

IT 0115-PL
One Needleless
Injection Port
Low Pressure
78in (198cm)



IV Extension Sets

Non-Pyrogenic, Sterile Fluid Path, Latex Free, DEHP Free


IT 0116-PL
Two Needle-less
Injection Ports
Low Pressure
78in (198cm)

IT 0122-PL
3-Way Stopcock, Needle-
less Injection Port
Low Pressure
78in (198cm)

IT 0606-PL
60 inch Extension
60in (152cm)

IT 0701-PL
3-Way Stopcock
Extension
10in (25cm)

IT 0360-PL
360in Coiled MRI
Extension Set
10in (25cm)



S1: Spike Assy w/open-close Vent

S1a: Protective Cap for Spike

S2: Spike w/open-close Vent

S2a: Protective Cap for Spike

C1: Roller Clamp

J1: Injection Port

L1: Female Luer Lock

L2: Male Luer Lock

C2: Pinch Clamp

J2: Needleless Injection Port

L1a: Protective Cap for Female Luer

L2a: Protective Cap for Male Luer

C3: Slide Clamp

J3: T-Injection Port

K1: 3-Way Stopcock

M2: Measured Volume Chamber

Features & Specifications

PL Medical IV Administration Sets and Extensions are fully customizable.

We can modify existing part numbers to accommodate your needs: such as adjusting lengths, adding or removing components etc.

Also we can create fully customized lines from the ground up.



PL Medical offers IV extension sets with exceptional quality and competitive value.

PL Medical I.V. extension sets are distributed in over 50 countries worldwide.

Our manufacturing facilities are both ISO 9001 and ISO 13485 certified.

Our I.V. extension sets have been manufactured with the invaluable inputs from scientists and research and development professionals throughout the world.

117 West Dudley Town Road
Bloomfield, CT 06002
info@plmedical.com
<http://www.plmedical.com>
T: (800) 874-0120

plmedical