

PL Medical offers IV Sets and Extension Lines with exceptional quality and competitive value.

PL Medical IV Sets and Lines are distributed in over 50 countries worldwide. Our manufacturing facilities are both ISO 9001 and ISO 13485 certified. Our IV Sets and Lines have been manufactured with the invaluable inputs from scientists and research and development professionals throughout the world.

I.V. EXTENSION SET – 30”


Part Number:
IT0109

Approximate priming volume – 4 ml

Fluid path sterile

Latex free

(1) Protective cover (2) Female adaptor (3) Tube (4) Roller clamp (5) Option lock (6) Male adaptor (7) Protective cap



I.V. EXTENSION SET – 30”


Part Number:
IT0110

Approximate priming volume – 4 ml

Fluid path sterile

Latex free

(1) Protective cover (2) Female adaptor (3) Tube (4) Slide clamp (5) Option lock (6) Male adaptor (7) Protective cap



I.V. EXTENSION SET – 42” 2 INJECTION SITE, MALE LUER LOCK ADAPTOR


Part Number:
IT0111

Approximate priming volume – 5.6 ml

Fluid path sterile

Latex free

(1) Cap (2) Female adaptor (3) Slide clamp (4) Tube (5) Injection site (6) Option lock (7) Male adaptor (8) Protective cap



I.V. EXTENSION SET – 5.25” WITH T AND INJECTION SITE


Part Number:
IT0112

Approximate priming volume – 0.7 ml

Fluid path sterile

Latex free

(1) Cap (2) Female adaptor (3) Tube (4) Slide clamp (5) T-connector (6) Protective cap (7) Injection site




INTERLINK SYSTEM – 6.5” T CONNECTOR I.V. EXTENSION SET

Part Number:
IT0113

Approximate priming volume – 0.6 ml
Fluid path sterile
Latex free

(1) Cap (2) Female adaptor (3) Tube (4) Slide clamp (5) T-connector (6) Protective cap (7) Injection site




I.V. EXTENSION LINE LOW PRESSURE 78” WITH ONE INJECTION SITE

Part Number:
IT0115

Approximate priming volume – 6 ml
Fluid path sterile
Latex free

(1) Protective cover (2) Female adaptor (3) Slide clamp (4) Cap (5) Injection site (6) Tube (7) Male adaptor (8) Protective cap




I.V. EXTENSION LINE 78” LOW PRESSURE WITH TWO INJECTION SITES

Part Number:
IT0116

Approximate priming volume – 6 ml
Fluid path sterile
Latex free

(1) Protective cover (2) Female adaptor (3) Slide clamp (4) Cap (5) Injection site (6) Tube (7) Male adaptor (8) Protective cap




I.V. EXTENSION LINE LOW PRESSURE 78” WITH 3 WAY STOPCOCK AND INJECTION SITE

Part Number:
IT0122

Approximate priming volume – 6 ml
 Fluid path sterile
 Latex free

(1) 3 way stopcock (2) Tube (3) Cap (4) Injection site (5) Male adapter (6) Protective cap (7) Cap




I.V. EXTENSION LINE LOW PRESSURE 78" WITH ONE INJECTION SITE

Part Number:
 IT0123

Approximate priming volume – 6 ml
 Fluid path sterile
 Latex free

(1) Protective cap (2) Female adapter (3) Tube (4) Slide clamp (5) Injection site (6) Male adapter (7) Protective cap (8) Cap



360" COILED MRI EXTENSION SET

Part Number:
 IT0360

Approximate priming volume – 12 ml
 Length 360"
 Fluid path sterile
 Fixed male and female luers
 Adjust tubing length by pulling from coiled tube as needed
 DEHP free
 Latex free

(1) Protective cap for male fitting adaptor (2) Male luer lock adaptor (3) DEHP free tubing (1.4 mm x 2.8 mm x 24") (4) "C" clamp (5) Female luer lock (6) Male luer lock adaptor (7) Coiled DEHP free tubing (1.4 mm x 2.8 mm x 243") (8) Female luer lock adaptor (9) Male luer lock adaptor (10) DEHP free tubing (1.4 mm x 2.8 mm x 93") (11) Female luer lock adaptor (12) Protective cap for male fitting adaptor

