

PL Medical offers surgical blades with exceptional value that puts them a cut above the rest.

PL Medical surgical blades are known for their precision. They assure reliability and durability to the surgeon. Our surgical blades have been manufactured with the invaluable inputs from scientists and research and development professionals all over the world.

- **Available in Sterile and Non-Sterile**
  - Guaranteed for 5 years—as long as the packaging remains unbroken
- **Available in Stainless Steel and Carbon Steel**
- **Box of 100**
- **Individually foil packed with V.C.I. paper**
- **Computer controlled grinding and polishing ensures a consistent sharp edge**
- **Heat treated to minimize breakage**



## Available Sizes

- #10 100Blades/Box, 100 Boxes/Case
- #11 100Blades/Box, 100 Boxes/Case
- #12 100Blades/Box, 100 Boxes/Case
- #12B 100Blades/Box, 100 Boxes/Case
- #15 100Blades/Box, 100 Boxes/Case
- #15C 100Blades/Box, 100 Boxes/Case
- #20 100Blades/Box, 100 Boxes/Case
- #21 100Blades/Box, 100 Boxes/Case
- #22 100Blades/Box, 100 Boxes/Case
- #23 100Blades/Box, 100 Boxes/Case
- #24 100Blades/Box, 100 Boxes/Case
- #25 100Blades/Box, 100 Boxes/Case
- #60 100Blades/Box, 100 Boxes/Case