



**JANUARY 31 - FEBRUARY 2, 2019
BOSTON CONVENTION & EXHIBITION CENTER- Booth 1551**

PL Medical is proud to be attending the Yankee Dental Congress 2019 in Boston. We will be located in Booth 1551. Please visit us to view our dosimetry service, scalpels, and blades.

Our Sales Manager, Martin Faber, will be in attendance to answer any questions you may have regarding the PL Medical dosimetry service and surgical blades. We look forward to meeting everyone.

About our Dosimetry Services:

PL Medical supplies accredited dosimetry services to industries in which radiation monitoring is desired or mandated by state and/or federal regulatory standards. PL Medical is accredited to provide dosimetry services by the National Voluntary Laboratory Accreditation Program (NVLAP Lab Code 200666-0), a division of National Institute of Standards and Technology (NIST).

- Reporting periods available: Monthly, Bi-Monthly, Quarterly & Semi-Annual
- Badge Description: 4 element dosimeters
- Minimum Reportable Dose: 10 mrem (0.10mSv)
- Energies Measured: X-Ray, Beta, Gamma
- Lower Limit of Detection: 6mrem (0.06 mSv)
- Useful Dose Range (10mrem – 1,000 rad)



Advantages:

- Customer Care – you will always be in contact with a “live person
- Environmentally Stable – Our TLD’s are rugged and generally hold-up to most day-to-day activities. If however you experience an issue where the badge is damaged we will replace the badge for you at no additional cost.
- 300mrem is automatically re-analyzed to ensure a correct reading
- Control the Additional Fees Others Charge: No account maintenance or set-up fees
- Flexibility of the Wear Period – We offer quarterly, monthly, bi-monthly and semi-annual exchanges. Pick a service that is right for your company. Additionally, you can start the service any time you wish, even in the middle of the month.
- Reports – Reports are provided to the customer 1-2 weeks after receipt of the badge at our laboratory. If you wish to have your report provided earlier, please contact us.
- Single point detector with chain of custody serial number – Each TLD is individually bar-coded and tracked from our laboratory.
- Each dose report is reviewed to check if it meets ALARA limits

About our Surgical Blades and Scalpels:

- Come in Stainless Steel or Carbon Steel
- Available in Sterile and Non-Sterile
- Have a variety of sizes 10,11,12,12B,15,15C,20,21,22,23
- Also have safety and thumb scalpels



Hours of Operations: M – F, 8:00 – 5:00 ET

Order Online: <http://www.plmedical.com>

Accredited by:



Panasonic ideas for life

A copy of our scope of accreditation can be found on <http://www.plmedical.com/index.php/accreditations/>